



INSTITUT KESIHATAN NEGARA
NATIONAL INSTITUTES OF HEALTH (NIH)
KEMENTERIAN KESIHATAN MALAYSIA
INSTITUT PENYELIDIKAN PERUBATAN (IMR)
Aras 6, Blok C6
No.1, Jalan Setia Murni U13/52
Seksyen U13 Setia Alam
40170 SHAH ALAM, SELANGOR

Telefoni : 603-3362 7900
Faksimili : 603-3362 7901

Ruj.kami: IMR/IDRC/VIRO(SA)/2301/23/01(131)

Tarikh: 3^{Mei} April 2021

YBrs. Dr,

LAPORAN UJIAN PENILAIAN KIT NOVEL CORONAVIRUS (COVID-19): NEWGENE BIOENGINEERING COVID-19 ANTIGEN DETECTION KIT

Dengan segala hormat merujuk kepada perkara di atas.

2. Adalah dimaklumkan bahawa pihak kami telah menjalankan ujian penilaian terhadap kit novel Coronavirus (2019-nCoV) *test kit* seperti yang dimohon oleh pihak YBrs. Dr.

3. Bersama-sama ini dikemukakan laporan bertajuk *Performance of NEWGENE Bioengineering COVID-19 Antigen Detection Kit, Manufacturer: New Gene (Hangzhou) Bioengineering Co., Ltd., China*, untuk perhatian dan rujukan pihak YBrs. Dr.

Sekian, terima kasih.

“PRIHATIN RAKYAT : DARURAT MEMERANGI COVID-19”

“BERKHIDMAT UNTUK NEGARA”

Saya yang menjalankan amanah,

(DR. HAJI TAHIR BIN ARIS)

MPM: 28418 NSR: 127511

Pengarah

Institut Penyelidikan Perubatan (IMR)

☎ 03 – 3362 8906

✉ tahir.a@moh.gov.my

S.k.:

En. Ahmad Shariff Bin Hambali

Ketua Eksekutif

Pihak Berkuasa Peranti Perubatan (*Medical Device Authority*)

Aras 5 & 6, Prima 9, Prima Avenue II, Blok 3547, Persiaran Apec, 63000 Cyberjaya, Selangor

**INSTITUTE FOR MEDICAL RESEARCH
KUALA LUMPUR**

Performance of Newgene Bioengineering COVID-19 Antigen Detection Kit

Intended use

This product is suitable for the qualitative detection of novel coronavirus in nasopharyngeal swab or sputum samples. It provides an aid in the diagnosis of infection with novel coronavirus.

Manufacturer

New Gene (Hangzhou) Bioengineering Co., Ltd., China.

Test Principle

This kit is an immunochromatographic membrane assay that uses highly sensitive monoclonal antibodies to detect nucleocapsid from protein from SARS-CoV-2. When the sample is added into the sample well, colloidal-gold conjugated with the monoclonal antibody against the nucleocapsid protein of SARS-CoV-2 absorbed in the reagent pad are dissolved and migrate along with the sample. If SARS-CoV-2 antigen is present in the sample, the complex of the anti-SARS-CoV-2 conjugate and the virus will be captured by the specific anti-SARS-CoV-2 monoclonal antibodies coated on the test line region (T). Absence of T line suggests a negative result. To serve as a procedural control, a red line will always appear in the control line region (C).

Test Kit

The evaluation was carried out using this kit with the lot number of 20210228-01, and the expiry date is Feb 27, 2022.

Instrument Used

NA

Reagent and Sample Preparation, Result Interpretation

Kindly refer to product package insert in the attachment

Sample Used

Sample type = Deep throat saliva (DTS)

Known Positive Sample= 50

Known Negative Sample= 50



Total samples used for analysis = 100

Performance Analysis

Test		Tested Kit Assay		Interpretation
		Positive	Negative	
IMR In-house Panel	SARS-CoV2 Positive (COVID-19)	48	2	Sensitivity= 96%
	SARS-CoV2 Negative (COVID-19)	0	50	Specificity = 100%

Comments

This kit was tested using DTS samples that had been tested positive using our RT-qPCR test system and the Ct values were ranged from 12 to 28. Saliva samples that had been tested negative using COVID-19 RT-qPCR test system and cultured samples of Influenza A H3, Influenza A pdm09, Adenovirus and Dengue were selected as negative panels.

Test performed by		Reviewed by	Approved by
			 3/5/22
Pn. E. Kavithambigai A/P Ellan Evaluator 1/Officer in charge	Pn. Hariyati Md. Ali Evaluator 2 /Officer in charge	Dr. Ravindran A/L Thayan Head of Virology Unit, IDRC, IMR	Dr. Hj. Tahir bin Aris Director of IMR

Disclaimer:

The content of this report does not imply that this product is endorsed or recommended by the Institute for Medical Research (IMR). This report may not be reviewed, abstracted, quoted, reproduced or translated, in part or in whole, without the prior written permission of IMR. IMR does not allow the use of this report or reference in any manner in the labeling or advertising of this product or kit package. However, the permission to import, export, or place in the market is subject to Medical Device Authority, Ministry of Health Malaysia's approval.